

Aldehyde dehydrogenase, mitochondrial

Catalog No: #AP79264



Package Size: #AP79264-1 50ug #AP79264-2 100ug #AP79264-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Aldehyde dehydrogenase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant Protein
Other Names	ALDM,ALDH class 2,ALDH-E2,ALDHI
Accession No.	P05091Gene name:ALDH2
Uniprot	P05091
GeneID	217;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only