

T-lymphocyte activation antigen CD80

Catalog No: #AP79265



Package Size: #AP79265-1 50ug #AP79265-2 100ug #AP79265-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-lymphocyte activation antigen CD80
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	CD28LG,CD28LG1,LAB7,Activation B7-1 antigen,BB1,CTLA-4 counter-receptor B7.1,B7,CD80
Accession No.	P33681Gene name:CD80
Uniprot	P33681
GenelD	941;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

(Microbial infection) Acts as a receptor for adenovirus subgroup B.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only