

Interleukin-22

Catalog No: #AP79267

Package Size: #AP79267-1 50ug #AP79267-2 100ug #AP79267-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-22
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	UNQ3099/PRO10096,ILTIF,ZCYTO18,Cytokine Zcyto18,IL-10-related T-cell-derived-inducible factor,IL-TIF
Accession No.	Q9GZX6Gene name:IL22
Uniprot	Q9GZX6
GeneID	50616;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Cytokine that contributes to the inflammatory response in vivo.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only