

Growth/differentiation factor 15

Catalog No: #AP79268

Package Size: #AP79268-1 50ug #AP79268-2 100ug #AP79268-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Growth/differentiation factor 15
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	30-300AA
Other Names	MIC1,PDF,PLAB,PTGFB,Macrophage inhibitory cytokine 1,NSAID-activated gene 1 protein,NSAID-regulated gene 1 protein,Placental TGF-beta,Placental bone morphogenetic protein,Prostate differentiation factor,MIC-1,NAG-1,NRG-1
Accession No.	Q99988Gene name:GDF15
Uniprot	Q99988
GeneID	9518;
Calculated MW	29.81
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Images



Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only