

## C-C chemokine receptor type 5

Catalog No: #AP79274



Package Size: #AP79274-1 50ug #AP79274-2 100ug #AP79274-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-C chemokine receptor type 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	CMKBR5,CHEMR13,HIV-1 fusion coreceptor,CD195
Accession No.	P51681Gene name:CCR5
Uniprot	P51681
GeneID	1234;
Tag Info	His
Formulation	50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

(Microbial infection) Acts as a receptor for human immunodeficiency virus-1/HIV-1.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only