

Galactoside-binding soluble lectin 13

Catalog No: #AP79287



Package Size: #AP79287-1 50ug #AP79287-2 100ug #AP79287-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galactoside-binding soluble lectin 13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	6-138AA
Other Names	PLAC8,Galectin-13,Placental tissue protein 13,Gal-13,PP13
Accession No.	Q9UHV8Gene name:LGALS13
Uniprot	Q9UHV8
GeneID	29124;
Calculated MW	14.63
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Binds beta-galactoside and lactose. Strong inducer of T-cell apoptosis.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only