

Endothelin-3

Catalog No: #AP79294

Package Size: #AP79294-1 50ug #AP79294-2 100ug #AP79294-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Endothelin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	118-238AA
Other Names	Preproendothelin-3,PPET3
Accession No.	P14138Gene name:EDN3
Uniprot	P14138
GeneID	1908;
Calculated MW	13.31
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Endothelins are endothelium-derived vasoconstrictor peptides.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only