

Fatty acid-binding protein, brain

Catalog No: #AP79298



Package Size: #AP79298-1 50ug #AP79298-2 100ug #AP79298-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Fatty acid-binding protein, brain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	Recombinant protein
Other Names	Blbp,Brain lipid-binding protein,Brain-type fatty acid-binding protein,Fatty acid-binding protein 7,BLBP,B-FABP
Accession No.	P51880Gene name:Fabp7
Uniprot	P51880
GeneID	12140;
Tag Info	His
Formulation	50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

B-FABP could be involved in the transport of a so far unknown hydrophobic ligand with potential morphogenic activity during CNS development. It is required for the establishment of the radial glial fiber system in developing brain, a system that is necessary for the migration of immature neurons to establish cortical layers.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only