

# Mesothelin

Catalog No: #AP79330



Package Size: #AP79330-1 50ug #AP79330-2 100ug #AP79330-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mesothelin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	Recombinant protein
Other Names	Erc,Mpf,Pre-pro-megakaryocyte-potentiating factor,Protein expressed in renal carcinoma
Accession No.	Q9ERA7Gene name:Msln
Uniprot	Q9ERA7
GeneID	60333;
Tag Info	His
Formulation	50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Megakaryocyte-potentiating factor (MPF) may potentiate megakaryocyte colony formation.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only