

Growth arrest and DNA damage-inducible protein GADD45 gamma

Catalog No: #AP79354

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP79354-1 50ug #AP79354-2 100ug #AP79354-3 1mg

Description

Product Name	Growth arrest and DNA damage-inducible protein GADD45 gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	CR6,DDIT2,Cytokine-responsive protein CR6,DNA damage-inducible transcript 2 protein,DDIT-2
Accession No.	O95257Gene name:GADD45G
Uniprot	O95257
GeneID	10912;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Involved in the regulation of growth and apoptosis. Mediates activation of stress-responsive MTK1/MEKK4 MAPKKK.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only