

GDNF family receptor alpha-1

Catalog No: #AP79355



Package Size: #AP79355-1 50ug #AP79355-2 100ug #AP79355-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GDNF family receptor alpha-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	25-429AA
Other Names	GDNFRA, RETL1, TRNR1, RET ligand 1, TGF-beta-related neurotrophic factor receptor 1
Accession No.	P56159 Gene name: GFRA1
Uniprot	P56159
GeneID	2674;
Calculated MW	44.55
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Receptor for GDNF. Mediates the GDNF-induced autophosphorylation and activation of the RET receptor.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only