

## Transmembrane glycoprotein NMB

Catalog No: #AP79356

Package Size: #AP79356-1 50ug #AP79356-2 100ug #AP79356-3 1mg


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Transmembrane glycoprotein NMB
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	322-486aa
Other Names	UNQ1725/PRO9925,HGFIN,NMB,Transmembrane glycoprotein HGFIN
Accession No.	Q14956Gene name:GPNMB
Uniprot	Q14956
GeneID	10457;
Calculated MW	18.15
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Images



## Background

Could be a melanogenic enzyme.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

---

Note: This product is for in vitro research use only