

Transmembrane glycoprotein NMB

Catalog No: #AP79357



Package Size: #AP79357-1 50ug #AP79357-2 100ug #AP79357-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transmembrane glycoprotein NMB
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	67-166AA
Other Names	Dchil,Hgfin,Nmb,DC-HIL,Dendritic cell-associated transmembrane protein,Osteoactivin
Accession No.	Q99P91Gene name:Gpnmb
Uniprot	Q99P91
GeneID	93695;
Calculated MW	11
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Could be a melanogenic enzyme.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only