

Protein jagged-2

Catalog No: #AP79358



Package Size: #AP79358-1 50ug #AP79358-2 100ug #AP79358-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein jagged-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	27-307aa
Other Names	null
Accession No.	Q9QYE5Gene name:Jag2
Uniprot	Q9QYE5
GeneID	16450;
Calculated MW	30.91
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Putative Notch ligand involved in the mediation of Notch signaling. Plays an essential role during limb, craniofacial and thymic development. May be involved in myogenesis and in the development of peripheral and central nervous systems.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only