

# Protein Wnt-4

Catalog No: #AP79361



Package Size: #AP79361-1 50ug #AP79361-2 100ug #AP79361-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein Wnt-4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	Recombinant protein
Other Names	Wnt-4
Accession No.	P22724Gene name:Wnt4
Uniprot	P22724
GenID	22417;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Ligand for members of the frizzled family of seven transmembrane receptors. May be an intracellular signaling molecule involved in segmentation of the forebrain. Is likely to signal over only few cell diameters (By similarity). Seems to be involved in kidney development.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only