

Testican-2

Catalog No: #AP79371



Package Size: #AP79371-1 50ug #AP79371-2 100ug #AP79371-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Testican-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	310-376AA
Other Names	UNQ269/PRO306,KIAA0275,TICN2,SPARC/osteonectin, CWCV, and Kazal-like domains proteoglycan 2
Accession No.	Q92563Gene name:SPOCK2
Uniprot	Q92563
GeneID	9806;
Calculated MW	7.37
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May participate in diverse steps of neurogenesis. Binds calcium.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only