

## DNA damage-regulated autophagy modulator protein

1

Catalog No: #AP79374

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP79374-1 50ug #AP79374-2 100ug #AP79374-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	DNA damage-regulated autophagy modulator protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	113-222AA
Other Names	DRAM,Damage-regulated autophagy modulator
Accession No.	Q8N682Gene name:DRAM1
Uniprot	Q8N682
GeneID	55332;
Calculated MW	12.1
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Lysosomal modulator of autophagy that plays a central role in p53/TP53-mediated apoptosis. Not involved in p73/TP73-mediated autophagy.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only