

DnaJ homolog subfamily C member 27

Catalog No: #AP79379

Package Size: #AP79379-1 50ug #AP79379-2 100ug #AP79379-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DnaJ homolog subfamily C member 27
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	122-231AA
Other Names	RABJS,RBJ,Rab and DnaJ domain-containing protein
Accession No.	Q9NZQ0Gene name:DNAJC27
Uniprot	Q9NZQ0
GeneID	51277;
Calculated MW	12.1
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

GTPase which can activate the MEK/ERK pathway and induce cell transformation when overexpressed. May act as a nuclear scaffold for MAPK1, probably by association with MAPK1 nuclear export signal leading to enhanced ERK1/ERK2 signaling.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only