

WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 2

Catalog No: #AP79380

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP79380-1 50ug #AP79380-2 100ug #AP79380-3 1mg

Description

Product Name	WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	Recombinant protein
Other Names	Gasp1,Growth and differentiation factor-associated serum protein 1,GASP-1
Accession No.	Q7TQN3Gene name:Wfikkn2
Uniprot	Q7TQN3
GeneID	278507;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Protease-inhibitor that contains multiple distinct protease inhibitor domains. Probably has serine protease- and metalloprotease-inhibitor activity (By similarity). Inhibits the biological activity of mature myostatin, but not activin.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only