

Adropin

Catalog No: #AP79382



Package Size: #AP79382-1 50ug #AP79382-2 100ug #AP79382-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Adropin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	34-76AA
Other Names	UNQ470/PRO830,C9orf165,Energy homeostasis-associated protein
Accession No.	Q6UWT2Gene name:ENHO
Uniprot	Q6UWT2
GeneID	375704;
Calculated MW	4.73
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Involved in the regulation of glucose homeostasis and lipid metabolism.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only