

Phospholipase D3

Catalog No: #AP79384



Package Size: #AP79384-1 50ug #AP79384-2 100ug #AP79384-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phospholipase D3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	260-490aa
Other Names	Choline phosphatase 3,HindIII K4L homolog,Hu-K4,Phosphatidylcholine-hydrolyzing phospholipase D3
Accession No.	Q8IV08Gene name:PLD3
Uniprot	Q8IV08
GeneID	23646;
Calculated MW	25.41
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May be involved in APP processing.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only