

## Serine/threonine-protein kinase 17B

Catalog No: #AP79390



Package Size: #AP79390-1 50ug #AP79390-2 100ug #AP79390-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein kinase 17B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	DRAK2,DAP kinase-related apoptosis-inducing protein kinase 2
Accession No.	O94768Gene name:STK17B
Uniprot	O94768
GeneID	9262;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Phosphorylates myosin light chains (By similarity). Acts as a positive regulator of apoptosis.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only