

## Beta-1,4-galactosyltransferase 6

Catalog No: #AP79400

Package Size: #AP79400-1 50ug #AP79400-2 100ug #AP79400-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Beta-1,4-galactosyltransferase 6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	Recombinant protein
Other Names	UDP-Gal:beta-GlcNAc beta-1,4-galactosyltransferase 6,UDP-galactose:beta-N-acetylglucosamine beta-1,4-galactosyltransferase 6
Accession No.	Q9WVK5Gene name:B4galt6
Uniprot	Q9WVK5
GeneID	56386;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Required for the biosynthesis of glycosphingolipids.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only