

CDYL2 Polyclonal Antibody

Catalog No: #40726

Package Size: #40726-1 50ul #40726-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	CDYL2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ICC
Species Reactivity	Hu
Specificity	CDYL2 Polyclonal Antibody detects endogenous levels of CDYL2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human CDYL2.
Target Name	CDYL2
Other Names	CDYL2; Chromodomain Y-like protein 2; CDY-like 2
Accession No.	Swiss-Prot: Q8N8U2NCBI Gene ID: 124359
Uniprot	Q8N8U2
GeneID	124359;
SDS-PAGE MW	56kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300. Immunocytochemistry: 1/200 - 1/1000.

Images



Western Blot analysis of HEPG2 cells using CDYL2 Polyclonal Antibody



Western Blot analysis of K562 Jurkat cells using CDYL2 Polyclonal Antibody

Note: This product is for in vitro research use only