

Utrophin

Catalog No: #AP79413



Package Size: #AP79413-1 50ug #AP79413-2 100ug #AP79413-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Utrophin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	DMDL,DRP1,Dystrophin-related protein 1,DRP-1
Accession No.	P46939Gene name:UTRN
Uniprot	P46939
GeneID	7402;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May play a role in anchoring the cytoskeleton to the plasma membrane.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only