

GDNF family receptor alpha-like

Catalog No: #AP79414



Package Size: #AP79414-1 50ug #AP79414-2 100ug #AP79414-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GDNF family receptor alpha-like
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	19-351aa
Other Names	UNQ9356/PRO34128,C6orf144
Accession No.	Q6UXV0Gene name:GFRAL
Uniprot	Q6UXV0
GeneID	389400;
Calculated MW	36.63
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only