

Multidrug resistance protein 1

Catalog No: #AP79415



Package Size: #AP79415-1 50ug #AP79415-2 100ug #AP79415-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Multidrug resistance protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	996-1280AA
Other Names	MDR1,PGY1,ATP-binding cassette sub-family B member 1,P-glycoprotein 1,CD243
Accession No.	P08183Gene name:ABCB1
Uniprot	P08183
GeneID	5243;
Calculated MW	31.35
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Energy-dependent efflux pump responsible for decreased drug accumulation in multidrug-resistant cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only