

DnaJ homolog subfamily B member 3

Catalog No: #AP79416



Package Size: #AP79416-1 50ug #AP79416-2 100ug #AP79416-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DnaJ homolog subfamily B member 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-145AA
Other Names	HCG3
Accession No.	Q8WWF6Gene name:DNAJB3
Uniprot	Q8WWF6
GeneID	414061;
Calculated MW	15.95
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

May operate as a co-chaperone of the male germ cell- and haploid stage-specific Hsp70 proteins.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only