

Cystatin-related protein 1

Catalog No: #AP79418

Package Size: #AP79418-1 50ug #AP79418-2 100ug #AP79418-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cystatin-related protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rat
Immunogen Description	Recombinant protein
Other Names	Crp1, Androgen-regulated 20 kDa protein, Prostatic 22 kDa glycoprotein P22K16/P22K20
Accession No.	P22282 Gene name: Andpro
Uniprot	P22282
GeneID	25030;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only