

CEP70 Polyclonal Antibody

Catalog No: #40728

Package Size: #40728-1 50ul #40728-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CEP70 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Mk
Specificity	CEP70 Polyclonal Antibody detects endogenous levels of CEP70 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CEP70.
Target Name	CEP70
Other Names	CEP70; BITE; Centrosomal protein of 70 kDa; Cep70; p10-binding protein
Accession No.	Swiss-Prot: Q8NHQ1NCBI Gene ID: 80321
Uniprot	Q8NHQ1
GeneID	80321;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using CEP70 Polyclonal Antibody

Note: This product is for in vitro research use only