

Cathelicidin-1

Catalog No: #AP79437



Package Size: #AP79437-1 50ug #AP79437-2 100ug #AP79437-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cathelicidin-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Immunogen Description	Recombinant protein
Other Names	BAC1,Bactenecin-1,Cyclic dodecapeptide,Bac1
Accession No.	P22226Gene name:CATHL1
Uniprot	P22226
GeneID	282164;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Potent microbicidal activity; active against S.aureus and E.coli.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only