

## UPF0556 protein C19orf10

Catalog No: #AP79444



Package Size: #AP79444-1 50ug #AP79444-2 100ug #AP79444-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	UPF0556 protein C19orf10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	32-173AA
Other Names	IL25,MYDGF,Interleukin-25,Stromal cell-derived growth factor SF20,IL-25
Accession No.	Q969H8Gene name:C19orf10
Uniprot	Q969H8
GeneID	56005;
Calculated MW	15.62
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only