

Granzyme K

Catalog No: #AP79449



Package Size: #AP79449-1 50ug #AP79449-2 100ug #AP79449-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Granzyme K
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	TRYP2,Fragmentin-3,Granzyme-3,NK-tryptase-2,NK-Tryp-2
Accession No.	P49863Gene name:GZMK
Uniprot	P49863
GeneID	3003;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

null

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only