

Tumor necrosis factor receptor superfamily member 5

Catalog No: #AP79450



Package Size: #AP79450-1 50ug #AP79450-2 100ug #AP79450-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tumor necrosis factor receptor superfamily member 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	20-193AA
Other Names	Tnfrsf5,B-cell surface antigen CD40,Bp50,CD40L receptor,CD40
Accession No.	P27512Gene name: Cd40
Uniprot	P27512
GeneID	21939;
Calculated MW	19.14
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Receptor for TNFSF5/CD40LG. Transduces TRAF6- and MAP3K8-mediated signals that activate ERK in macrophages and B cells, leading to induction of immunoglobulin secretion.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only