

## Cytokine receptor common subunit gamma

Catalog No: #AP79459



Package Size: #AP79459-1 50ug #AP79459-2 100ug #AP79459-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytokine receptor common subunit gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	Interleukin-2 receptor subunit gamma,gammaC,p64,IL-2 receptor subunit gamma,CD132
Accession No.	P31785Gene name:IL2RG
Uniprot	P31785
GeneID	3561;
Tag Info	His
Formulation	50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Common subunit for the receptors for a variety of interleukins.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only