

Antibacterial protein PR-39

Catalog No: #AP79469



Package Size: #AP79469-1 50ug #AP79469-2 100ug #AP79469-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Antibacterial protein PR-39
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Immunogen Description	30-130AA
Other Names	null
Accession No.	P80054Gene name:PR39
Uniprot	P80054
GeneID	407610;
Calculated MW	11.11
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Exerts a potent antimicrobial activity against both E.coli and B.megaterium.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only