

## Kallikrein 1-related peptidase b22

Catalog No: #AP79470



Package Size: #AP79470-1 50ug #AP79470-2 100ug #AP79470-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kallikrein 1-related peptidase b22
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Mouse
Immunogen Description	CYEDKYNIWLGKNKLFQDEPSAQ
Other Names	Klk-22, Klk22, Beta-NGF-endorpeptidase, Epidermal growth factor-binding protein type A, Glandular kallikrein K22, Nerve growth factor beta chain endopeptidase, Tissue kallikrein 22, EGF-BP A, mGK-22
Accession No.	P15948 Gene name: Klk1b22
Uniprot	P15948
GeneID	13646;
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Glandular kallikreins cleave Met-Lys and Arg-Ser bonds in kininogen to release Lys-bradykinin.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only