

Hyaluronan mediated motility receptor

Catalog No: #AP79474



Package Size: #AP79474-1 50ug #AP79474-2 100ug #AP79474-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hyaluronan mediated motility receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	180-312AA
Other Names	IHABP,RHAMM,Intracellular hyaluronic acid-binding protein,Receptor for hyaluronan-mediated motility,CD168
Accession No.	O75330Gene name:HMMR
Uniprot	O75330
GeneID	3161;
Calculated MW	14.63
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Involved in cell motility. When hyaluronan binds to HMMR, the phosphorylation of a number of proteins, including PTK2/FAK1 occurs. May also be involved in cellular transformation and metastasis formation, and in regulating extracellular-regulated kinase (ERK) activity.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only