

Hepatocyte growth factor receptor

Catalog No: #AP79475



Package Size: #AP79475-1 50ug #AP79475-2 100ug #AP79475-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hepatocyte growth factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	HGF/SF receptor,Proto-oncogene c-Met,Scatter factor receptor,Tyrosine-protein kinase Met,SF receptor
Accession No.	P08581Gene name:MET
Uniprot	P08581
GeneID	4233;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Acts as a receptor for Listeria internalin inIB, mediating entry of the pathogen into cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only