

Argininosuccinate synthase

Catalog No: #AP79482



Package Size: #AP79482-1 50ug #AP79482-2 100ug #AP79482-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Argininosuccinate synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1-240AA
Other Names	ASS, Citrulline--aspartate ligase
Accession No.	P00966 Gene name: ASS1
Uniprot	P00966
GeneID	445;
Calculated MW	26.4
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.) Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Is indirectly involved in the control of blood pressure.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only