

Protein-arginine deiminase type-2

Catalog No: #AP79487



Package Size: #AP79487-1 50ug #AP79487-2 100ug #AP79487-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein-arginine deiminase type-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	438-665AA
Other Names	KIAA0994,PAD2,PDI2,PAD-H19,Peptidylarginine deiminase II,Protein-arginine deiminase type II
Accession No.	Q9Y2J8Gene name:PADI2
Uniprot	Q9Y2J8
GeneID	11240;
Calculated MW	25.08
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Catalyzes the deimination of arginine residues of proteins.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only