

Interleukin-17D

Catalog No: #AP79491



Package Size: #AP79491-1 50ug #AP79491-2 100ug #AP79491-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-17D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	Recombinant protein
Other Names	UNQ3096/PRO21175,Interleukin-27,IL-27
Accession No.	Q8TAD2Gene name:IL17D
Uniprot	Q8TAD2
GeneID	53342;
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Induces expression of IL6, CXCL8/IL8, and CSF2/GM-CSF from endothelial cells.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only