

CGI-43 Polyclonal Antibody

Catalog No: #40735

Package Size: #40735-1 50ul #40735-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CGI-43 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Mk
Specificity	CGI-43 Polyclonal Antibody detects endogenous levels of CGI-43 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human CGI-43.
Target Name	CGI-43
Other Names	CCZ1; C7orf28A; CGI-43; Vacuolar fusion protein CCZ1 homolog; CCZ1B; C7orf28B; Vacuolar fusion protein CCZ1 homolog B; Vacuolar fusion protein CCZ1 homolog-like
Accession No.	Swiss-Prot: P86791/P86790NCBI Gene ID: 221960
Uniprot	P86790
GeneID	221960;51622;
SDS-PAGE MW	56kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using CGI-43 Polyclonal Antibody

Note: This product is for in vitro research use only