

# Secretory carrier-associated membrane protein 1

Catalog No: #AP79495



Package Size: #AP79495-1 50ug #AP79495-2 100ug #AP79495-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Secretory carrier-associated membrane protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	2-155AA
Other Names	SCAMP
Accession No.	O15126Gene name:SCAMP1
Uniprot	O15126
GeneID	9522;
Calculated MW	16.94
Tag Info	His
Formulation	50mM NaH <sub>2</sub> PO <sub>4</sub> , 500mM NaCl Buffer with 500mM Imidazole, 10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

## Background

Functions in post-Golgi recycling pathways. Acts as a recycling carrier to the cell surface.

## References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only