

Death-inducer obliterator 1

Catalog No: #AP79499

Package Size: #AP79499-1 50ug #AP79499-2 100ug #AP79499-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Death-inducer obliterator 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Human
Immunogen Description	1725-2034AA
Other Names	C20orf158,DATF1,KIAA0333,Death-associated transcription factor 1,DATF-1
Accession No.	Q9BTC0Gene name:DIDO1
Uniprot	Q9BTC0
GeneID	11083;
Calculated MW	34.1
Tag Info	His
Formulation	50mM NaH ₂ PO ₄ , 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)
Storage	Store at -20C. (Avoid repeated freezing and thawing.)Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Background

Isoform 2: Displaces isoform 4 at the onset of differentiation, required for repression of stemness genes.

References

Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.

Note: This product is for in vitro research use only