

## Ferrous fumarate

Catalog No: #S0001


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ferrous fumarate
Brief Description	Inhibitors
Purification	90.00%
Target Name	Others
Calculated MW	169.9
Formulation	C <sub>4</sub> H <sub>2</sub> FeO <sub>4</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Metabolism system  
 SMILES: C(=C/C(=O)[O-])C(=O)[O-].[Fe+2]  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: >280°C  
 Solubility: Soluble in Water  
 Appearance: Brownness Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.  
 Stock solution can be stored below -20 °C for several month.

## References

1. Garca-Casal MN, Leets I. J Food Sci. 2014 Apr;79(4):H706-12.

Note: This product is for in vitro research use only