

Etamsylate

Catalog No: #S0002


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Etamsylate
Brief Description	Inhibitors
Purification	99.55%
Target Name	Thrombin inhibitor
Calculated MW	263.31
Formulation	C10H17NO5S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
SMILES: OC1=CC(=C(O)C=C1)S(O)(=O)=O
Pathways: Others
Receptor: Thrombopoiesis
Boiling pt: Melting pt: 125°C
Solubility: Water: 53 mg/mL (201.28 mM); DMSO: 53 mg/mL (201.28 mM)
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Garay RP, et al. Am J Ther, 2006, 13(3), 236-247.

Note: This product is for in vitro research use only