

Methylene Blue trihydrate

Catalog No: #S0006


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Methylene Blue trihydrate
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others inhibitor
Calculated MW	373.9
Formulation	C16H24ClN3O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Other SMILES:O.O.O.[Cl-].CN(C)C1=CC2=C(C=C1)N=C1C=CC(C=C1S2)=[N+](C)C Pathways:Others Receptor:Others Boiling pt:Melting pt:190C?Solubility:Water: 43,600 mg/L.Appearance:Brownness Powder Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Rager T, et al. Phys Chem Chem Phys. 2012 Jun 14;14(22):8074-82.

Note: This product is for in vitro research use only