

Atenolol

Catalog No: #S0007


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Atenolol
Brief Description	Inhibitors
Purification	99.89%
Target Name	Adrenergic Receptor agonist
Calculated MW	266.34
Formulation	C14H22N2O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular System
SMILES: CC(C)NCC(O)COC1=CC=C(CC(N)=O)C=C1
Pathways: GPCR/G Protein Receptor: α -adrenergic receptor
Boiling pt: 508C? Melting pt: 154C
Solubility: soluble in Water
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Nuttall SL, et al. J Clin Pharm Ther. 2003 Jun;28(3):179-86.

Note: This product is for in vitro research use only