

## Atenolol

Catalog No: #S0007


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Atenolol
Brief Description	Inhibitors
Purification	99.89%
Target Name	Adrenergic Receptor agonist
Calculated MW	266.34
Formulation	C14H22N2O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area: Cardiovascular System  
SMILES: CC(C)NCC(O)COC1=CC=C(CC(N)=O)C=C1  
Pathways: GPCR/G Protein Receptor:  $\alpha$ -adrenergic receptor  
Boiling pt: 508°C  
Melting pt: 154°C  
Solubility: soluble in Water  
Appearance: White Powder  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Nuttall SL, et al. J Clin Pharm Ther. 2003 Jun;28(3):179-86.

Note: This product is for in vitro research use only